


Simple and straightforward Data analysis



Carbonaceous Aerosol Analysis Tool - CAAT

KEY FEATURES:

- **IMPORT** your data from Magee Scientific instrument
- **VALIDATE** your data.
- **ANALYZE** your validated data by using several procedures: Diurnal profile, Time Average, Absorption Coefficient, Angstrom exponent, Attenuation, Source Apportionment.
- **DISPLAY** your data by using various chart options.



Developed and Manufactured in
Europe by Aerosol d.o.o., Slovenia

Product specifications

Analyze your Magee Scientific instrument data with CAAT Software.

BENEFITS

Be assured that your data is correct and consistent

- Use CAAT as your primary tool for your analysis of data obtained from Magee Scientific instruments.
- Immediately check the proper “data status”.


Start using CAAT Software tool

- No additional requirements.
- Straightforward data evaluation and display.
- User friendly controls.





EASY INSTALLATION

Login or Register at [Magee Scientific User's group](#):

1. Search for CAAT product.
2. Go to SOFTWARE tab and download installation files.
3. Search for the downloaded files in your computer. Unzip to a selected destination folder.
4. Click on the “CAAT-Setup.exe” and follow the guided procedure.
5. You can start using the CAAT application immediately after installation is completed.



Scan the code for more info



GENERAL INQUIRIES:

Aerosol d.o.o., Kamniška 41, SI-1000 Ljubljana, Slovenia
tel: +386(1)4391700 fax: +386 59 191 221 www.aerosol.eu
US, Mexico, Canada:

Magee Scientific Corp., 1916 M.L. King Jr. Way,
Berkeley CA 94704, USA tel: +15108452801 fax: +15108457137
www.mageescientific.com

or the distributor responsible for your country.